

APPLANATION TONOMETER



Glass Agencies are manufacturer of Applanation Tonometer Under EROSE Brand.

□ SPECIFICATIONS:-

Measuring Range	0 to 10.64kpa (0 mmHg to 80mmHg)
Movement of light Circle	1.53*2=3.06mm
Diameter of Probe	7 MM
Moving Range of Probe	3 MM
Weight	300 Gms
Size	190x80x80mm
Manufacturer	Glass Agencies
Model	EROSE ATSL 0.0001
Brand	EROSE
Material	Metal
Packing	Box
Colour	Black
Usage	Clinic, Hospital, Medical Purpose

*Designs and details can be changes without prior notice for the purposes of improvement.



GLASS AGENCIES

5309, Anaj Mandi, Ambala Cantt. 133001 (Haryana) INDIA

Office : +91 89304-07858

Mobile : +91 98968-07858

E-mail : glassagencies@gmail.com

Website : www.glassagencies.com