

EROSE ASPHERIC LENS 20 D



Glass Agencies are manufacturer of
DOUBLE ASPHERIC LENS 20 D Under
EROSE Brand.

Specification :

Usage : Clinic, Hospital

Power/Resolution : 20 D

Material : Aluminium Ring, Glass lens

Manufacturer Name : Glass Agencies

Brand : EROSE

Color : Black

Shape : Round

Application : Laser Lens

Packaging Type : Box

This 20D in every physicians pocket. Our improved design delivers superior image clarity throughout the viewing field.

Computer optimized lens surfaces eliminate optical aberrations. A true 20 diopter design. Available in colors : Black, Blue, purple, Red Or Green

GLASS AGENCIES

5309, Anaj Mandi, Ambala Cantt. 133001 (Haryana) INDIA

Office : +91 171 2633027 Mobile : +91 98968-07858

E-mail : glassagencies@gmail.com

Website : www.glassagencies.com

