

EROSE ASPHERIC LENS 78 D



Glass Agencies are manufacturer of ASPHERIC LENS 78 D Under EROSE Brand.

Specification :

Usage : Clinic, Hospital

Power/Resolution : 78 D

Material : Aluminium Ring, Glass lens

Manufacturer Name : Glass Agencies

Brand : EROSE

Color: Black, Blue, Purple, Red or Green

Field of View : 84 Degree

Shape : Round

Application : Laser Lens

Packaging Type : Box

This Slit lamp lens offers an ideal combination magnification and field of view. It also providing an extended working distance from cornea. Superior glass option offers high clarity images. Ideal general diagnostic and treatment lens.



GLASS AGENCIES

5309, Anaj Mandi, Ambala Cantt. 133001 (Haryana) INDIA

Office : +91 171 2633027 Mobile : +91 98968-07858

E-mail : glassagencies@gmail.com

Website : www.glassagencies.com