

EROSE ASPHERIC LENS 90 D



Glass Agencies are manufacturer of ASPHERIC LENS 90 D Under EROSE Brand.

Specification :

Usage: Clinic, Hospital

Power/Resolution: 90 D

Material: Aluminium Ring, Glass lens

Brand: EROSE

Color : Black, Blue, Purple, Red or Green

Shape: Round

Application: Laser Lens

Field of View : 94 Degree

Packaging Type: Box

The standard slit lamp biomicroscopy lens. Our 90D design deliver superior clarity, stereopsis and outstanding small pupil capabilities. Compare it to your 90D and See the difference.

This Slit Lamp lens is the most popular power for non-contact fundus examination. It offers a wide-field fundus imaging out to and beyond the equator. A high magnification provides an increased illuminated field area which makes it ideal for posterior pole scanning.



GLASS AGENCIES

5309, Anaj Mandi, Ambala Cantt. 133001 (Haryana) INDIA

Office : +91 171 2633027 Mobile : +91 98968-07858

E-mail : glassagencies@gmail.com

Website : www.glassagencies.com