

Matrix Ophthalmics, S. Korea



Inspiring Ophthalmology....

# MATRIX AUTOREFRACTOMETER : R<sup>20+</sup> & MATRIX AUTOREFKERATOMETER:KR<sup>20+</sup>



**GLASS AGENCIES**  
5309, Anaj Mandi,  
Ambala Cantt. 133001  
(Haryana) INDIA  
Office : +91 171 2633027  
Mobile : +91 98968-07858

## SALIENT FEATURE:-

- ☛ New generation hi-resolution Image processor with 21st century ARM processor and Japanese optical path technology, resulting sharpest image with fastest measurement.
- ☛ Dual CCD camera technology with 5.7 inch hi-resolution color LCD display
- ☛ Accurate and wide measuring range with Auto & Manual mode of measurement (selectable).
- ☛ Provides both corneal curvature in mm and corneal power in Diopter ( available in KR20+)
- ☛ Scientifically designed Scenic Fixation target with auto fogging technology can even measure smallest pupil of 2mm size.
- ☛ On-Screen display of 10 measurements per eye with PD & average data for subjective test.
- ☛ Latest SMPS power supply system, allows working perfect even at low voltage or fluctuating supply voltage.
- ☛ Compatible with MATRIX CV 7000 auto vision tester (Auto Phoropter)