

# BINOCULAR INDIRECT OPHTHALMOSCOPE, FOUR FILTER

Glass Agencies are manufacturer of BINOCULAR INDIRECT OPHTHALMOSCOPE Under EROSE Brand.

## SPECIFICATION :

Usage : Clinic, Hospital, Laboratory

Automatic Grade : Manual

Color : Black

Filters : Cobalt Blue, Green ,Red and White

Manufacturer Name : Glass Agencies

Brand : EROSE

Handling : Portable

Material : Fiber

Type : Rechargeable Indirect Ophthalmoscope

- Has All Pupil Feature
- Brilliant, Natural & Even illumination from a 6v 20w Halogen Lamp.
- Stereo Optical System
- Built in Blue ,Green, Red and White filters & three step field of illumination.
- Compact and light weight for longer periods without fatigue.
- Precision optics and elegant mechanism, synchronizes the convergence and parallel adjustments.
- Also available in Fiber Optic Light system
- Interpupillary Distance adjustable between 50-75mm
- Head Band Indirect Ophthalmoscope
- Simple Control for Adjusting tshe Head Band

## STANDARD ACCESSORIES

- Carrying Case
- Co- Observation Teaching Mirror
- S. Depressor
- Fundus Chart



## GLASS AGENCIES

5309, Anaj Mandi, Nr B. D. High School, Ambala Cantt. 133001 (Haryana) INDIA

Office : +91 171 2633027 Mobile : +91 98968-07858

E-mail : [glassagencies@gmail.com](mailto:glassagencies@gmail.com)

Website : [www.glassagencies.com](http://www.glassagencies.com)

