

EROSE BINOCULAR COAXIAL MICROSCOPE



STANDARD ACCESSORIES

- 10X or 15x Eyepieces
- Dust Cover
- Empty Box for objective
- Dust Cloth

Glass Agencies are manufacturer
of **BINOCULAR COAXIAL**
MICROSCOPE Under **EROSE** Brand.

Specification :

Usage : Laboratory, Hospital , Clinic
Binocular Head : 360° Rotatable & 45° In
Magnification : 50X to 1000X / 75 X to 150
Manufacturer Name : Glass Agencies
Model : EROSE BCM 0.000 I
Material : Metal & Plastic
Colour : White & Black,
Operation Mode : Manual
LED Light : 3 Watts with Low / High Dimmer
Eye pieces : 10x / 15x E Wide Field (18 mm F
Objective : 5x, 10x, 40x, 100x (Oil)
Packaging Type : Thermocoal Box
Focusing System : Coaxial
Mechanical Stage : Double layered with knobs
Filter Holder : Swing out Filter Holder
Knob : for separate Coarse & Fine Adjustment
Condenser : Movable double lens Bright Field abb
condenser with iris diaphragm

GLASS AGENCIES

5309, Anaj Mandi, Ambala Cantt – 133001, Haryana, INDIA

Mobile : +91 9416024836, 98968-07858

E-mail : glassagencies@gmail.com

Website : www.glassagencies.com

