

EROSE LAB BINOCULAR MICROSCOPE

Glass Agencies are manufacturer
of **LAB BINOCULAR MICROSCOPE**
Under **EROSE** Brand.



Specification :

Usage : Laboratory, Hospital , Clinic
Binocular Head : 360° Rotatable & 45° In
Magnification : 50X to 1500X
Manufacturer Name : Glass Agencies
Model : EROSE LBM 001
Material : Metal & Plastic
Colour : White & Black,
Operation Mode : Manual
LED Lamp : 3 Watts with Low / High Dimmer
Eye pieces : 10x & 15x Wide Field
Objective : 5x, 10x, 40x, 100x (Oil)
Packaging Type : Thermocoal Box
Mechanical Stage : Movable
Filter Holder : Swing out Filter Holder
Knob : for separate Coarse & Fine Adjustment
Condenser : Movable double lens Bright Field abt
condenser with iris diaphragm

STANDARD ACCESSORIES

- 10X & 15x Eyepieces
- Dust Cover
- Empty Box for objective
- Dust Cloth

GLASS AGENCIES

5309, Anaj Mandi, Ambala Cantt – 133001, Haryana, INDIA

Mobile : +91 9416024836, 98968-07858

E-mail : glassagencies@gmail.com

Website : www.glassagencies.com

