

# EROSE DOCTOR MODEL OPHTHALMIC CHAIR UNIT



Glass Agencies are  
manufacturer of  
Doctor Model  
Ophthalmic Chair Unit

## Features

- This Model for optical shop and for the small rooms.
- Instrument smooth rotation and well-balanced top.
- Table top for Auto Refractometer and for the Lensometer for easy refraction
- Soft touch safety switch to prevent accidental chair movement.
- Over head task light that is bright, Focused and rheostat-Controlled.
- Centrally integrated control panel with touchstone keys.
- Movements across all four directions are operated silently by CE certified German DC motor.
- Sturdy pole to sustain projector, phoropter and Slit Lamp.
- Excellent maneuverability to offer best height and position of near vision holder and projector arm.
- In this chair, only chair going up and down, reclining is manual from chair back cover
- Only back position is reclining and foot rest not come up.
- Multifaceted controls from up/down movement & memory of zero position.
- Offers in built connectivity for near vision holder and projector.
- Highly flexible chair that rotates and reclines with manually.
- Foot switch to facilitate smooth control of upward and downward movement.
- Polyurethane seats that are anti static and comfortable for the patient.
- Electro mechanical movements that enable all four movements
- Heavy load-uplifting capacity with minimum power consumption.
- Fully integrated Headrest
- BRAND : EROSE
- Manufacturer : Glass Agencies

## GLASS AGENCIES

5309, Anaj Mandi, Ambala Cantt. 133001 (Haryana) INDIA

Office : +91 89304-07858

Mobile : +91 98968-07858

E-mail : [glassagencies@gmail.com](mailto:glassagencies@gmail.com)

Website : [www.glassagencies.com](http://www.glassagencies.com)

