

FIBRE OPTIC HALOGEN LIGHT SOURCE



➤ Glass Agencies are manufacturer of Fibre Optic Halogen Light Source Under EROSE Brand.

❖ Features:

- Dual Light System, 2 Lamp of 250 Watts
- Bulb Changer Switch
- Colour Temperature approx. 3400K
- In-built Heat Filter
- Manual Light Intensity Adjustment
- Manufacturer Name GLASS AGENCIES
- Double Fan Cooling System
- Metal Body Material
- In built Voltage Control Rectifier
- Inbuilt Safety Fuse
- Tungsten Halogen Type with 250K Reflector
- EROSE Brand
- Power supply 230V (AC) / 50Hz ± 10%
- Weight 6 Kg

GLASS AGENCIES

5309, Anaj Mandi, Ambala Cantt. 133001
(Haryana) INDIA

Office : +91 89304-07858

Mobile : +91 98968-07858

E-mail : glassagencies@gmail.com

Website : www.glassagencies.com

