

# FOUR MIRROR GONIOSCOPE



Glass Agencies are manufacturer of **FOUR MIRROR GONIOSCOPE** Under **EROSE** Brand.

Specification :

Type : Four Mirror Gonioscope

Usage : Hospital, Clinic

Shape : Round

Packaging Type : Box

Material : Glass

Brand : EROSE

Portable : Yes

Mirror Angle : 62 Degree

Frame Material : Aluminium

Manufacturer Name : Glass Agencies



Four Mirror angled at 62 Degree help in viewing all Quadrants of the anterior chamber without Rotation. The small contact Surface of 9mm diameter makes the lens very convenient to use it has a unique detachable handle. Perfect for fast and effective anterior segment examinations.

**GLASS AGENCIES**

5309, Anaj Mandi, Ambala Cantt. 133001 (Haryana) INDIA

Office : +91 171 2633027 Mobile : +91 98968-07858

E-mail : [glassagencies@gmail.com](mailto:glassagencies@gmail.com)

Website : [www.glassagencies.com](http://www.glassagencies.com)

