

FUME HOOD



➤ Glass Agencies are manufacturer of Fume Hood Under EROSE Brand.

Fume exhaust hoods are designed to contain and exhaust toxic, obnoxious Or otherwise harmful gases, vapors, mist etc to protect personal and equipment. A motor driven blower crates negative pressure within the chamber extracting the Contaminated air from work-area and expel it in to the atmosphere. The hood is ventilated of the air drawn from inside the laboratory which is used effectively for thorough sweep- out of work chamber.

The standard features 'EROSE' Fume Hood are:-

- Inbuilt centrifugal Exhaust unit
- Fluorescent tubes with safety electrical
- Internally located electric plugs
- Fibre Glass Sink Best Quality (Small size)
- Gas lines Standard, vacuum & Air connections (optional)
- Counter balanced sliding up & down Front Sash door
- Water tap & waste pipe outlet.

Utility chamber or storage unit beneath the fume hood unit to keep utility material such as glassware, material handling equipment and accessories with storage capacity starting From approx. 400 liter up to 1200 ltr. Approx. optional.

GLASS AGENCIES

5309, Anaj Mandi, Ambala Cantt. 133001

(Haryana) INDIA

Office : +91 89304-07858

Mobile : +91 98968-07858

E-mail : glassagencies@gmail.com

Website : www.glassagencies.com



FUME HOOD



Stainless steel chamber:- Fabricated of stainless steel sheet from inside and outside. (S.S-304) angular & pipe frame work, double chamber With strong centrifugal blower exhaust system.

Electro-galvanized steel chamber:- Fabricated of thick pre-electro-galvanized steel, havi duty sheet 20 or 18 gauge frame work, angular & pipe double chamber, Thoroughly epoxy powder-coated from inside and outside. With strong centrifugal exhaust system

Wooden Chamber:- Fabricated out of best quality thick weather resistant duro board clad with super white mica lamination from outside, inner chamber finished in chemical-resistant Epoxy paint and work surface covered with thick non-glossing & non rating Stainless steel or ceramic ties.

To work on 220/230 V AC 50 HZ

NOTE:-The inner chamber can be lined with polypropylene or Teflon Led sheet (optional) Following are the standard dimensions available-

Stainless Steel	Mild Steel	Wooden	Size of Working Chamber
EROSE FHSS-01	EROSE FHMS-01	EROSE FHW-01	900mm x 750mm x 1200mm
EROSE FHSS-02	EROSE FHMS-02	EROSE FHW-02	1200mm x 750mm x 1200mm
EROSE FHSS-03	EROSE FHMS-03	EROSE FHW-03	1500mm x 750mm x 1200mm

Other sizes can be designed as per customer specifications.
External Exhaust with various capacities available.(Cost Extra)

Note: All the Optional Accessories cost extra.

GLASS AGENCIES

5309, Anaj Mandi, Ambala Cantt. 133001
(Haryana) INDIA

Office : +91 89304-07858

Mobile : +91 98968-07858

E-mail : glassagencies@gmail.com

Website : www.glassagencies.com

