

ILLUMINATED EROSE TRIAL CASE WITH WOODEN BOX

Glass Agencies are manufacturer of ILLUMINATED TRIAL CASE WITH WOODEN BOX Under EROSE Brand.

Specification :

Application: USED FOR EYE TESTING

Usage: Hospital, Clinic

Types: Full Diameter Trial Set

Manufacturer Name : Glass Agencies

Brand : EROSE

Colour: Golden & Silver

Material: Metal and Glass

Number Of Lenses: 225

Packaging Type: Box

- Reduced aperture lenses fitted in metal mounts with Powder Coated.
- Deluxe Model
- +/- sph 32 pairs each.
- +/- cyl 20 pairs each.
- 10 Prisms, 7 Accessories.
- Trial Frame



GLASS AGENCIES

5309, AnajMandi, AmbalaCantt. 133001 (Haryana) INDIA

Office : +91 171 2633027

Mobile : +91 98968-07858

Email : glassagencies@gmail.com

Website : www.glassagencies.com