

EROSE MANUAL KERATOMETER



Glass Agencies are manufacturer of **Manual Keratometer** under **EROSE** Brand.

Specification:

- One position instrument.
- Measures both meridians without changing optical system.
- Positive fixation, Permits rapid measurement of central corneal area.
- Precision objectives, Assures durability, Achromatic lenses.
- Dual eye-level sighting system. Facilitates horizontal alignment.
- Two-way adjustable head and chin rests.
- Simple vertical adjustment to fit patient.
- Coincidence focusing system.
- Positive measurements to .02mm accuracy, 15x wide-angle piece.
- Astigmatic axis 0 to 180 Degree
- The achromatic optical systems to precise measurement of the corneal radius curvature and accurate measurement of radius curvature of contact lens
- Measuring range : radius of curvature.9.4-6.4 mm (in 0.05 mm steps)
- Diopter 36-52 (in 0.25 steps)
- Contents: 2 spare Bulb, Duster Cover
- Manufacturer Name : Glass Agencies
- Brand : EROSE
- Usage : Hospital, Clinic
- Type : Ophthalmic Optical Equipments
- Power Source : Electric
- Color : Black
- Packaging Type : Box

GLASS AGENCIES

5309, Anaj Mandi, Ambala Cantt. 133001 (Haryana) INDIA

Office : +91 89304-07858

Mobile : +91 98968-07858

E-mail : glassagencies@gmail.com

Website : www.glassagencies.com

