

LENSOMETER DOUBLE TARGET CROSS



Glass Agencies are manufacturer of
**LENSOMETER DOUBLE TARGET
CROSS** Under **EROSE** Brand.

Specification :

Usage: Hospital, Laboratory

Color: White

Material: Metal and Plastic

Brand: EROSE

Power Source : Electric

Portable: Yes

Operation Mode: Manual

Packaging Type: Box

- Crosslines with corona target
- External reading with tilling facility
- +/- to 25 Dioptors power reading
- 0 to 5 prism reading
- 0° to 180° cylindrical axis
- Battery and mains operated

GLASS AGENCIES

5309, Anaj Mandi, Ambala Cantt. 133001 (Haryana) INDIA

Office : +91 171 2633027

Mobile : +91 98968-07858

E-mail : glassagencies@gmail.com

Website : www.glassagencies.com