

# MILK CENTRIFUGE



➤ Glass Agencies are  
Manufacturer of Milk  
Centrifuge Under EROSE  
Brand.

## **SPECIFICATION:-**

- \* Electricity Operated Fat Testing Centrifuge Machine
- \* Inclusive of rotor to accommodate 8 / 12 butyrometers.
- \* Timer (0-5 Minutes)
- \* Brand Name EROSE
- \* Manufacturer Name Glass Agencies
- \* Mechanical Brake
- \* Anti-vibration rubber pad.
- \* Max Speed 1400 RPM

**GLASS AGENCIES**

5309, Anaj Mandi, Ambala Cantt. 133001 (Haryana) INDIA

Office : +91 89304-07858

Mobile : +91 98968-07858

E-mail : [glassagencies@gmail.com](mailto:glassagencies@gmail.com)

Website : [www.glassagencies.com](http://www.glassagencies.com)

