

# EROSE OBJECTIVE OIL IMMERSION 100 X LENS



Glass Agencies are  
manufacturer of **Objective  
Oil Immersion 100 X  
Lens** under **EROSE Brand**.

## Specification:

- Material : Brass Or Plastic
- Magnification : 100X
- Usage/Application : Microscope
- Manufacturer Name : Glass Agencies
- Brand : EROSE
- Shape : Round
- Colour : Black & Silver
- Portable : Portable

## Features:

- This is a 100X achromatic oil-immersion objective lens for DIN-standard microscopes.
- Intended primarily for biological applications, it is designed to work with microscopes using 160mm tube-length and 45mm parfocal distance.
- Achromatic correction provides improved color accuracy and resolution.

## GLASS AGENCIES

5309, Anaj Mandi, Ambala Cantt. 133001 (Haryana) INDIA

Office : +91 89304-07858

Mobile : +91 98968-07858

E-mail : [glassagencies@gmail.com](mailto:glassagencies@gmail.com)

Website : [www.glassagencies.com](http://www.glassagencies.com)

