

SINGLE MIRROR GONIOSCOPE

Glass agencies are manufacturer of Single Mirror Gonioscope under EROSE Brand.

Due to its small size a Single Mirror Gonioscope is very convenient to handle. The mirror angle is 62 degree which when rotated through 360 degree gives a clear view of the anterior chamber angle. Its small contact surface makes it very useful while diagnosing adults as well as children with small palpebral fissure.



* Made in glass

* Contact medium required

SPECIFICATION :

* Single Mirror

* MIRROR ANGLE : 62 degree

* MANUFACTURER : Glass Agencies

* BRAND : EROSE



GLASS AGENCIES

5309, Anaj Mandi, Ambala Cantt.133001 (Haryana) INDIA

Off. : +91 171 2633027, Mob. : +91 9896807858

Email : glassagencies@gmail.com

Website : www.glassagencies.com