

SYNOPTHOPHORE



SPECIFICATION :

Glass Agencies are manufacturer of SYNOPTTHOPHORE Under EROSE Brand.

Movement of Optical Tubes

HORIZONTAL

Adduction + 50°

Adduction - 40°

VERTICAL

Hyper 30°

Hypo 30°

TORSIONAL

Incylo 20°

Excyclo 20°

Slide Illumination : Rheostat controlled 12V Lamp for each Slide. After Image Illumination by Lamp (for better illumination).

Brand : EROSE

Manufacturer Name : Glass Agencies

Auto Flashing : Auto Flashing of Slide illumination either SIMULTANEOUS or ALTERNATE in RAPID & VARIABLE Modes.

Mode & Mode Selection : Five Modes of Slide illumination namely NORMAL, FLASHING RIGHT, FLASHING LEFT, FLASHING R+L & AUTO FLASHING, can be selected by a single selector knob.

Optional Accessories : Haidinger Brushes. A Set of Slides Containing 9 Pairs, Motorized Table.



GLASS AGENCIES

5309, Anaj Mandi, Ambala Cantt. 133001 (Haryana) INDIA

Office : +91 171 2633027 Mobile : +91 98968-07858

E-mail : glassagencies@gmail.com

Website : www.glassagencies.com