

TEAT OPENER



➤ Glass Agencies are manufacturer of Teat Opener Under EROSE Brand.

➤ SPECIFICATION:-

- Material Stainless Steel
- Use In Animal Treatment
- Usage Hospital
- Length 118 mm
- Colour Silver
- Product Type Standard Veterinary Instrument
- Feature Comfortable to use & Good in Quality
- Properties High Strength, Good Quality
- Brand EROSE
- Manufacturer Glass Agencies
- Teat Opener is used to open a Collide Canal.

GLASS AGENCIES

5309, Anaj Mandi, Ambala Cantt. 133001 (Haryana)
INDIA

Office : +91 89304-07858

Mobile : +91 98968-07858

E-mail : glassagencies@gmail.com

Website : www.glassagencies.com

