

TEAT MILK SYPHON



□ Glass Agencies are manufacturer of Teat Milk Siphon Under EROSE Brand.

❖ SPECIFICATION:-

▪Material	Stainless Steel
▪Use In	Animal Treatment
▪Usage	Hospital
▪Colour	Silver
▪Length	120 mm
▪Product Type	Standard Veterinary Instrument
▪Feature	Comfortable to use & Good in Quality
▪Properties	High Strength, Good Quality
▪Brand	EROSE
▪Manufacturer	Glass Agencies

A Milk **Siphon** is a specialized instrument used in veterinary medicine and dairy farming for taking milk samples directly from the **teat** canal of animals.

GLASS AGENCIES

5309, Anaj Mandi, Ambala Cantt. 133001 (Haryana)
INDIA

Office : +91 89304-07858

Mobile : +91 98968-07858

E-mail : glassagencies@gmail.com

Website : www.glassagencies.com

