

# THREE MIRROR GONIOSCOPE



Glass Agencies are manufacturer of **THREE MIRROR GONIOSCOPE** Under **EROSE** Brand.

## Specifications :

Usage : Hospital, Clinic

Shape : Round

Manufacturer Name : Glass Agencies

Brand : EROSE

Material : Glass

Mirror Angle : 59, 66, 73 Degree

This lens has three mirror angle at 59 Degree, 66 Degree and 73 Degree and arranged at 120 Degree intervals to permit viewing of the fundus and anterior chamber angle. The posterior pole is also observable through the central axis of the lens.

## GLASS AGENCIES

5309, Anaj Mandi, Ambala Cantt. 133001 (Haryana) INDIA

Office : +91 171 2633027 Mobile : +91 98968-07858

E-mail : [glassagencies@gmail.com](mailto:glassagencies@gmail.com)

Website : [www.glassagencies.com](http://www.glassagencies.com)

