

Two Mirror Gonioscope

Glass agencies are manufacturer of Two Mirror Gonioscope under EROSE BRAND.

Two Mirror Gonioscope is designed for viewing the anterior chamber angle and central retina. Two opposing mirrors inclined at 62 degree provide a complete 360 degree view of angle with only 90 degree rotation of the lens. Similar to a single mirror gonioscope, the small contact surface helps in examination of children and adults with small fissures.

- * Made in glass
- * Contact Medium Required

SPECIFICATION :

- * Two Mirror
- * Mirror Angle : 62/62 degree
- * Manufacturer : Glass Agencies
- * BRAND : EROSE



GLASS AGENCIES

5309, Anaj Mandi, Ambala cantt- 133001 (Haryana) INDIA

Off. : +91 1712633027 , Mob. : +91 9896807858

Email : glassagencies@gmail.com,

Website : www.glassagencies.com

