

CONVENTIONAL MANUAL ROTARY MICROTOME

- Glass Agencies are manufacturer of **Conventional Manual Rotary Microtome** under EROSE Brand.

High Precise Manual Rotary Microtome Is Wildly Used In Clinical Pathology Lab And Educational Purpose

Model	EROSE CMRM 0.0001
Operation	Manual
Section Thickness Range	1µm to 50µm in Increment of 1 µm)
Knife Length	120mm
Knife Length (optional)	140mm/160mm
Specimen Orientation	8° x,y axis 360°arround z axis
Vertical Stroke	55mm
Horizontal Advancement	20mm
Maximum Sectioning Area	32mm x 25mm
Disposable Blade Holder	optional
Dimension (DXWXH)	300mm x 300mm x 290mm (approx)
Weight	30 kg (Approx.)

- Accurate Specimen Advancement
- Safety Mechanism For Specimen Advancement.
- Smooth Hand Wheel With Locking System.
- Compatible With Disposable Blade Attachment.

GLASS AGENCIES

5309, Anaj Mandi, Ambala Cantt. 133001 (Haryana) INDIA

Office : +91 89304-07858

Mobile : +91 98968-07858

E-mail : glassagencies@gmail.com

Website : www.glassagencies.com

