

# FULLY AUTOMATIC CRYOSTAT



- Glass Agencies are manufacturer of **Fully Automatic Cryostat** under EROSE Brand.

## ❖ Specification:

- Automatic sectioning with motorized feed drive.
- Fast freezing cooling shelf  $-45^{\circ}\text{C}$ .
- Dual compressor touch control display.
- Manufacturer Glass Agencies
- Brand EROSE
- Section thickness  $1\mu\text{m}-100\mu\text{m}$ .
- Adjustable sectioning speed & retraction Disposable holder with anti roll plate.
- Weekly programmable start & stop time.

**GLASS AGENCIES**

5309, Anaj Mandi, Ambala Cantt. 133001 (Haryana) INDIA

Office : +91 89304-07858

Mobile : +91 98968-07858

E-mail : [glassagencies@gmail.com](mailto:glassagencies@gmail.com)

Website : [www.glassagencies.com](http://www.glassagencies.com)

