

# GROSSING WORK STATION

- Glass Agencies are manufacturer of **Grossing Work Station** under EROSE Brand.

Grossing Work Station With Back Draft Ventilation Made Up Of High Grade 30 Stainless Steel This High End Instrument Is Equid With Continuous Air Flow Over The Entire Work Are To Remove The Formalin Fumes.

## Stranded Accessory:

- Inbuilt & continues Ventilation System
- Magnetic Tool Holder
- Document Holder
- Hot & cold Water Faucet
- Dissecting Board

## Optional Accessory:

- Globus Holder
- Towel Holder
- Sturdy Perforated Work Table
- Magnifying Glass
- Monitor & Camera Imaging System
- Sprinkle System With Adjustable Flow
- Spot Light
- Integrated Illumination System

## GLASS AGENCIES

5309, Anaj Mandi, Ambala Cantt. 133001 (Haryana) INDIA

Office : +91 89304-07858

Mobile : +91 98968-07858

E-mail : [glassagencies@gmail.com](mailto:glassagencies@gmail.com)

Website : [www.glassagencies.com](http://www.glassagencies.com)

