

VACUUM TISSUE PROCESSOR

- Glass Agencies are manufacturer of **Vacuum Tissue Processor** under **EROSE Brand**.

Microprocessor Controlled Carousel Tissue Processor Provide Gentle And Highly Reliable Specimen Processing. For Faster And Substantially Improved Specimen Processing, This High End Instrument is Equipped With Vacuum Unit For All Wax And Reagent Station.

Technical Specifications

Model	EROSE VTP 0.002
Operation	Automated with Microprocessor controlled
S.S Basket	1 No
Capacity	100 Tissue Cassettes
Reagent Beaker	10 Nos.
Paraffin Bath	2 Nos, Electronically controlled
Beaker capacity	2 Ltr
Paraffin bath capacity	2 Ltr
Safety device	Yes
Delay function	Yes
Vacuum condition	Yes / Selected
Programme storage	9 Programme
Agitation	Rotary as well as Up/Down
Weight	28 Kg (Approx)

Instrument can be supplied with 9 reagent beaker and 3 paraffin bath Station.

- Broad Digital Display.
- Feather Touch Key Pad.
- Battery Back Up To Retains Data.
- Cycle Time Can Be Feed From 1 Minute To 999 Minute For Each Station.
- Second Specimen Transport Basket Can Be Load To Process Double Of Tissue Sample.
- 2 Or 3 Electronically Controlled Paraffin Bath,.
- Selectable Vacuum Condition For All 12 Station.

GLASS AGENCIES

5309, Anaj Mandi, Ambala Cantt. 133001 (Haryana) INDIA

Office : +91 89304-07858

Mobile : +91 98968-07858

E-mail : glassagencies@gmail.com

Website : www.glassagencies.com

