

# BROCK STRING



➤ Glass Agencies are manufacturer of Brock String Under EROSE Brand.

## ■ SPECIFICATION:-

❖ Brand	EROSE
❖ Rop Colour	White
❖ Handle Material	Plastic
❖ Length	12 Feet
❖ Manufacturer	Glass Agencies
❖ Beads Material	Plastic
❖ Number of beads	3
❖ Model Number	EROSE OBS 0.001

A Brock String is a common vision therapy tool. It consists of a white flexible cord or string, 12 feet in length, with colored beads that can be moved to various positions along the length of the string. The Brock offers instant feedback to the participants if their eyes are working together to focus on an object at various distances. The Brock string is often used in vision therapy for patients with vergence disorders such as Convergence Insufficiency, patients with Amblyopia (also called Lazy Eye) and for patients with Strabismus.

## GLASS AGENCIES

5309, Anaj Mandi, Ambala Cantt. 133001 (Haryana) INDIA

Office : +91 89304-07858

Mobile : +91 98968-07858

E-mail : [glassagencies@gmail.com](mailto:glassagencies@gmail.com)

Website : [www.glassagencies.com](http://www.glassagencies.com)

