


# EROSE DIGITAL BLOOD CELL COUNTER



Glass Agencies are manufacturer of **Digital Blood Cell Counter** under EROSE Brand.

## ❖ Specification :

- Display: : 16 x 2 LINES ALPHANUMERIC LCD.
- Key Pad : Total 16 Keys 4 Function & 16 Category Keys
- Each Category : 000 to 999 Counts
- Total Counts : 0000 to 9999
- Power Supply : 230 V + 10 %, 50 Hz.
- Battery : 12 V Battery (Optional)
- Manufacturer Name : Glass Agencies
- BRAND : EROSE
- Mode : Count Mode , Finished Mode, Percentage Mode
- The instrument consists of Twelve keys. The function of some keys is as follows:-
  - 1) ALL CLEAR : This key is used to start fresh counting & Clear all data.
  - 2) % : This key is used to view percentage.
  - 3) DEC : This key is used to view total counts.
  - 4) ENT : This key is used to remove last entry made by mistake.
- EROSE Blood Cell Counter incorporate latest microcontroller based technology to pack features of differential counting and percentage calculation in single instrument to ease blood cell counting procedure.

## GLASS AGENCIES

5309, Anaj Mandi, Ambala Cantt. 133001 (Haryana) INDIA

Office : +91 89304-07858

Mobile : +91 98968-07858

E-mail : [glassagencies@gmail.com](mailto:glassagencies@gmail.com)

Website : [www.glassagencies.com](http://www.glassagencies.com)

