


DIGITAL LABORATORY CENTRIFUGE



➤ Glass Agencies are manufacturer of Digital Laboratory Centrifuge Under EROSE Brand.

■ SPECIFICATION:-

❖ Brand	EROSE
❖ Type	Benchtop
❖ Display	Digital Speed Indicator, Digital Timer
❖ Manufacturer	Glass Agencies
❖ Min / Max Speed	500 / 5000 RPM
❖ Model Number	EROSE DLC 0.001
❖ Rotor Size	Angle head , Test Tube - 6, 8 (Optional 12, 16 & 20 Test Tubes)
❖ Stepless Sped Regulator	

EROSE DLC 0.001 is suitable for routine sample analysis in Medical, Hospital & Industrial Laboratories. It have digital Speed indicator & 0-99 Minute digital countdown timer. It has safety lid interlock to prevent cover opening during centrifugation. Heavy Duty CRC construction finished with powder coating.

GLASS AGENCIES

5309, Anaj Mandi, Ambala Cantt. 133001 (Haryana) INDIA

Office : +91 89304-07858

Mobile : +91 98968-07858

E-mail : glassagencies@gmail.com

Website : www.glassagencies.com

