

DIGITAL ULTRASONIC CLEANER



➤ Glass Agencies are manufacturer of Digital Ultrasonic Cleaner Under EROSE Brand.

SPECIFICATION:-

Tank Volume	1.8L	2L	2.5L	3.5L	4L	5.5L
Inner Dimension	150×137×110	150×135×145	235×135×100	295×145×100	235×135×100	300×153×150
Ultrasonic Wattage	50W	50W	100W	100W	100W	150W
Tank Volume	6L	9L	12L	16L	20L	25L
Inner Dimension	300×240×100	300×240×150	500×300×100	300×300×200	500×300×150	500×300×200
Ultrasonic Wattage	200W	200W	300W	300W	300W	600W

Display LED display with timer
 Digital Timer 5 To 60 minutes
 Inner / Outer Tank Material 304 Grade Stainless Steel

- *Uniform cavitations through out the tank.
- *Using 40 Khz Frequency for effective cleaning with low noise.
- *In a very short time it removes all dirt & foreign bodies from components which are immersed in the cleaning liquid.

GLASS AGENCIES

5309, Anaj Mandi, Ambala Cantt. 133001
 (Haryana) INDIA

Office : +91 89304-07858

Mobile : +91 98968-07858

E-mail : glassagencies@gmail.com

Website : www.glassagencies.com

