

EROSE VITRECTOMY LENS SET

Glass Agencies are
Manufacturer of
**VITRECTOMY
LENS SET** under
EROSE Brand.



❖ SPECIFICATIONS:

• This vitrectomy lens set comprises of a set of 7 contact lenses and 1 ring to hold the contact lenses in position. The main advantage of this system is that it minimizes the dependence of a well-trained assistant.

- Manufacturer Name : Glass Agencies
- Brand : EROSE
- High quality glass optics for clear images.
- Unique silicon storage cum sterilization container.
- Autoclavable lenses.
- Suturable metal ring.
- Flexibility to order individual lenses and rings.

GLASS AGENCIES

5309, Anaj Mandi, Ambala Cantt. 133001 (Haryana) INDIA

Office : +91 89304-07858

Mobile : +91 98968-07858

E-mail : glassagencies@gmail.com

Website : www.glassagencies.com

