

VISCOSITY / FORD CUP WITH STAND



➤ Glass Agencies are manufacturer of Viscosity Cup With Stand Under EROSE Brand.

➤ SPECIFICATION:-

- | | |
|--|---|
| ▪ Usage | Industrial and Laboratory |
| ▪ Material | Cup is made from Bras, Stand is made from Cast Iron |
| ▪ Viscosity | B1, B2, B3,B4 & B6 |
| ▪ Weight | 800 Gms |
| ▪ Volume | 100ml |
| ▪ Brand | EROSE |
| ▪ Manufacturer | Glass Agencies |
| ▪ It Suitable for measuring kinematic viscosity. Flow time of liquids. | |
| ▪ The cups can be supplied separately or with an adjustable stand | |

GLASS AGENCIES

5309, Anaj Mandi, Ambala Cantt. 133001
(Haryana) INDIA

Office : +91 89304-07858

Mobile : +91 98968-07858

E-mail : glassagencies@gmail.com

Website : www.glassagencies.com

