

EROSE LABORATORY BUNSEN BURNERS



Glass Agencies are manufacturer of **Laboratory Bunsen Burners** under EROSE Brand.

- **Specification :**
 - Usage : Laboratory
 - Material : Cast Iron & Brass
 - Color : Any
 - Base Size : 80mm
 - Manufacturer Name : Glass Agencies
 - Brand : EROSE
 - Shape : Round
-
- Cast iron base powder coated.
 - Brass pipe chrome plated.
 - Stop cock to control Gas Pressure.
 - Pipe Dia 12 mm with brass air controller to adjust flame.
 - Bunsen Burner with Stopcock for use with natural gas, offers durability, precision, and safety for laboratory heating applications. Its compact design, stopcock control, and compatibility with natural gas make it an essential tool for laboratories seeking a versatile and dependable heat source for scientific experiments and procedures.

GLASS AGENCIES

5309, Anaj Mandi, Ambala Cantt. 133001 (Haryana) INDIA

Office : +91 89304-07858

Mobile : +91 98968-07858

E-mail : glassagencies@gmail.com

Website : www.glassagencies.com

