

# TABLET DISINTEGRATION



➤ Glass Agencies are manufacturer of Tablet Disintegration Apparatus Under EROSE Brand.

Disintegration Tester is a solid state instrument designed for the accurate estimation of disintegration time of tablets as per IP/USP standards. The instrument is designed to test One / Two batches of six tablets. simultaneously. The unit is extremely useful for pharmaceutical industries. A motorized shaft moves in up and down direction at a frequency of 30 cycles per minute, holds 6 glass tubes of standard size 7.5 mm x 20 mm with stainless steel wire gauge bottom, which dip in water in a beaker of 1000 ml. Cylinder moves up and down with standard dish inside in a liter beaker filled of distilled water which is maintained at 37 C by thermostatically controlled hot plates (without beaker). The instrument has One / two basket assemblies with independent timers and motors. Each timer can be set and used individually. The basket speed is accurately controlled and maintained as per IP, USP specifications.

- A Single Test Apparatus.
- B Double Test Apparatus

- Optional
  - Digital Temp Display Meter
  - Digital R.P.M Display Meter.

## GLASS AGENCIES

5309, Anaj Mandi, Ambala Cantt. 133001  
(Haryana) INDIA

Office : +91 89304-07858

Mobile : +91 98968-07858

E-mail : [glassagencies@gmail.com](mailto:glassagencies@gmail.com)

Website : [www.glassagencies.com](http://www.glassagencies.com)

