

# TISSUE EMBEDDING STATION



- Glass Agencies are manufacturer of **Tissue Embedding Station** under EROSE Brand.

## ❖ Specification:

- Microprocessor controlled five independent heating zone.
- Setting Temp. Range: 50-70°C
- Automatic system start and stop (7day memory)
- Manufacturer Glass Agencies
- Brand EROSE
- Temperature display panel for all heated areas.
- Volume of the paraffin: 3 to 5ltrs.
- Two tissue preservation boxes, paraffin nozzle, optional cold spot.

**GLASS AGENCIES**

5309, Anaj Mandi, Ambala Cantt. 133001 (Haryana) INDIA

Office : +91 89304-07858

Mobile : +91 98968-07858

E-mail : [glassagencies@gmail.com](mailto:glassagencies@gmail.com)

Website : [www.glassagencies.com](http://www.glassagencies.com)

