

SEMI AUTOMATIC MICROTOME

- Glass Agencies are manufacturer of **Semi Automatic Microtome** under **EROSE Brand**.

Semi-Automatic Microtome Has A High Precession Stepper Motor Driven And Microprocessor Controlled Spaceman Advancement System For Uniform Advancement And Precise Positioning Of Spaceman

Model	EROSE SAM0.0003
Operation	Semi Automatic
Section Thickness Range	0.5µm To 100µm
Increment Value	From 0.5µm To 5µm in 0.5µm From 5µm To 20µm in 1µm From 20µm To 60µm in 5µm From 60µm To 100µm in 10µm
Trimming Thickness Range	0 To 500µm (with The Increment Of 5µm 10µm And 50µm)
Spaceman Retraction Setting Range	Sum To 100µm (With The Increment Of 5µm)
Specimen Orientation	8XY Axis 360 around 2 Axis
Vertical Stroke	60mm
Horizontal Advancement	25mm
Maximum Sectioning Area	50mm X 50mm
Dimension (DXWXH)	460mm x 520mm x 320mm (approx)
Weight	24 kg (Approx.)
Remote Control.	Optional

Instrument Will Be Supplied With Standard Object Clamp.

High Precession Motorizes Spaceman Advancement System

Lcd Display With Feather Touch Key Pad.

Counter Weight Balanced Hand Wheel And Ergonomically Located Hand Wheel Lock

Acoustic Or Acroscopic Signal Of Specimen Feed At End And Home Position..

.Low And High Profile Compatible Blade Holder

Blade Holder With integrated Finger Protection Guard And Blade Ejector

Spacious Waste Tray With Arm Rest.

GLASS AGENCIES

5309, Anaj Mandi, Ambala Cantt. 133001 (Haryana) INDIA

Office : +91 89304-07858

Mobile : +91 98968-07858

E-mail : glassagencies@gmail.com

Website : www.glassagencies.com

