



Binocular Head having Hard Coated Prisms  
Bright Field horizontal illuminator with centering provision  
Quadruple revolving nose piece  
Special Graduated mechanical stage  
Variable light control through Heavy-Duty transformer

**Objectives : 10x, 45x**  
**Optional objectives : 5x, 100x**  
**Eye-pieces : WF10x (paired)**  
**Magnification : 100x to 450x.**

**Packaging:** Inner Thermocol/Styrofoam Packing and  
Outer Cardboard Box.

# Binocular Metallurgical Microscope



## GLASS AGENCIES

5309, ANAJ MANDI, NR B. D. HIGH SCHOOL, AMBALA CANTT. 133001 (HARYANA) INDIA

MOBILE : +91 98968-07858, +91 94160-24836

E-MAIL : GLASSAGENCIES@GMAIL.COM,

WEBSITE : WWW.GLASSAGENCIES.COM